

Development of a myocardial perfusion phantom

Gijs de Vries, s1854526

Revision 0.100

Detailed Design

2

Development of a myocardial perfusion phantom

3

G.J. de Vries, s1854526

4

Tuesday 26th February, 2019

5 **Preface**

6 [todo] this

7 G.J. (Gijs) de Vries

8 Enschede, 13th of February 2019

9 **Contents**

10	1 Introduction	1
11	2 Research methodology	2
12	2.1 Main research question	2
13	2.2 Concept of Operations	2
14	2.3 Requirements and Architecture	2
15	2.4 Detailed Design	3
16	3 Concept design	4
17	3.1 Global restraints	4
18	3.2 Design categories	4
19	3.3 Concepts	7
20	3.4 Comparison	11
21	3.5 Definitive choice	13
22	4 Detailed design	14
23	A Appendix: Mind map	15
24	Bibliography	16

1 Introduction

Myocardial Perfusion Imaging (MPI), or, simply put, the imaging of the blood flow in the heart muscle, plays an important role in diagnosing heart failure or detecting Coronary Artery Disease (CAD). Imaging systems like Computed Tomography (CT), Magnetic Resonance Imaging (MRI), Single-Photon Emission Computed Tomography (SPECT), or Positron Emission Tomography (PET) can visualise a (radioactive) contrast bolus in the supplying arteries and in underlying myocardial tissue, whose flow can give an indication of narrowed or blocked blood vessels.

Many variations in the visualisation process of myocardial perfusion, including variations in hard- and software, can (significantly) influence the outcome and in turn have consequences for patient treatment. These variations need to be validated against a well-known baseline.

A myocardial perfusion phantom will be developed that is able to simulate the blood flow in the heart muscle, i.e. the myocardium, and is able to mimic cardiac defects like (significant) stenosis.

Document overview

[todo] This section

Abbreviations

CAD Coronary Artery Disease

CT Computed Tomography

HLA Horizontal Longitudinal Axis

MPI Myocardial Perfusion Imaging

MRI Magnetic Resonance Imaging

PET Positron Emission Tomography

SA Short Axis

SPECT Single-Photon Emission Computed
Tomography

VLA Vertical Longitudinal Axis

2 Research methodology

This chapter serves as a summary of the previously answered research questions and gives an overview of the research questions to come.

2.1 Main research question

Can patient treatment reliably depend on the D-SPECT, using dynamic scanning, in myocardial perfusion imaging?

Answer

As of Tuesday 26th February, 2019, the main research question has not been answered.

Answered in

It will be answered in the final report of the master's thesis.

Based on

The answer will be based on the developed myocardial perfusion phantom and the experiments performed at the ZGT Hengelo.

2.2 Concept of Operations

Is the D-SPECT's dynamic scanning, in comparison with other modalities(CT, MRI, PET, or SPECT), suitable for quantitative perfusion imaging?

What must the myocardial perfusion phantom be able to simulate?

Answer

The D-SPECT is relatively new in the Netherlands, but it is more widely employed in Japan, Canada, France, and Great-Britain. The highly specialised nature (for cardiac purposes), the patient friendly design, the ability to scan faster and more accurate at significant dose reductions, make the D-SPECT suitable for quantitative myocardial perfusion.

The myocardial perfusion phantom will have to simulate a patient, with stenotic artery (or arteries), in a physiological way, which is compatible with clinical protocol and software.

Answered in

The answer to this research question can be found in the system requirements document section 2.2.1 and 2.2.2, respectively.

Based on

The answer to this research question is based on the literature review and background research performed in the project plan, chapter 2.

2.3 Requirements and Architecture

What are the requirements for a myocardial perfusion phantom that can be used in combination with commonly used clinical software?

Answer

The requirements are specified in tables corresponding in the system requirements document.

87 **Answered in**

88 The answer to this research question can be found in the system requirements document,
89 chapters 2 and 3.

90 **Based on**

91 The answer to this research question is based on interviews with a part of the direct stakehold-
92 ers, as specified in the project plan section 3.2.

93 **2.4 Detailed Design**

94 *How can the myocardial perfusion phantom meet the clinical requirements and mimic the per-*
95 *fusion of a human heart?*

96 **Answer**

97 It will be answered in the this detailed design document.

98 **Answered in**

99 The answer can be found in chapter 4.

100 **Based on**

101 The answer to this research question is based on a mind map which results in different concept.
102 The most promising concept, based on the requirements, is developed further into a detailed
103 design.

3 Concept design

This chapter defines the global restraints on the concepts and presents an overview of the different concepts.

3.1 Global restraints

This section describes the global restraints on the concepts, as stated in the system requirements.

3.1.1 Myocardium

To ensure compatibility to the clinical software, *4DM*, the myocardium's HLA and VLA must have the shape of a horseshoe and the SA must have the shape of a circle. *4DM* requires these shapes to determine the contours and consequently determine the myocardial flow.

3.1.2 Modality

The main research question is based around the, relatively new in the Netherlands, D-SPECT's dynamic scanning. Therefore, the modality is bounded to the D-SPECT. As mentioned in section 2.2, the D-SPECT is a suitable choice for myocardial perfusion imaging but still requires validation, which is the goal of the PhD research of which this project is part of.

3.1.3 Tracer

The tracer and injection method are fixed due to clinical, *4DM*, and dynamic scanning requirements. The clinical (and D-SPECT) protocol use ^{99m}Tc (Technetium) Tetrofosmin for myocardial perfusion imaging. To ensure proper dynamic scanning results, the tracer must be injected with a pump. The pump can repeatedly inject tracer with identical volumes and injection speeds.

3.1.4 Flow Type

Based on background information on the *4DM* software, a decision has been made to use a non-pulsatile flow. The D-SPECT uses gated measurements such that images are extracted at the same point in time of the cardiac cycle. Furthermore, initial experiments with the D-SPECT, with non-pulsatile flow, have been performed on February 5, 2019. The prototype set-up from the individual project, with a dialysis tube, was placed in the D-SPECT and the TAC was extracted. This curve showed proper similarity to TACs extracted from patients.

3.2 Design categories

Each concept is based on the mind map shown in appendix A. The following sections describe each aspect of the five main categories, respectively:

- Myocardium,
- Modality,
- Design,
- Tracer, and
- Flow.

3.2.1 Myocardium

The myocardium has two subcategories, the shape and the model.

Shape

As described in section 3.1.1, the shape of the myocardium is especially important for the 4DM software since it searches for the contours of the left ventricle's walls. Therefore, the shape of the myocardium is fixed (as per system requirements):

- Vertical Longitudinal Axis (VLA): Horseshoe (figure).
- Horizontal Longitudinal Axis (HLA): Horseshoe.
- Short Axis (SA): circle.




Figure 3.1: Left ventricle's myocardial shapes: VLA, HLA, SA, respectively (Niwaz, 2015)

Model

[todo] this section, depending on 4dm.

The 4DM software offers different segment models in addition to the full 17-segment model, i.e. The goal is to develop a phantom with 17 segments. However, 4DM's ability to downscale to a fewer-segment model, provides the opportunity for a more simplified first prototype.

3.2.2 Modality

As described in section 3.1.2, the modality is fixed to the D-SPECT's dynamic scanning.

3.2.3 Design

The phantom can be designed in three different ways: using a 1-, 2-, or 4-chamber design. The 1-chamber design simulates only the left ventricle with the myocardium. The 2-chamber design simulates either the left and right ventricles, or the left atrium and ventricle. The other combinations, i.e. left ventricle and right atrium, and any combination without the left ventricle, do not hold any additional benefits. The 4-chamber design contains all heart chambers to simulate the flow as physiological as possible. A 3-chamber design is not considered since it has no physiological structure nor does it have an added benefit over a 1- or 2-chamber design.

3.2.4 Tracer

As described in section 3.1.3, the tracer protocol is fixed.

3.2.5 Flow

Generator

The flow in the phantom can be realised by a various of methods (or combinations thereof):

- Peristaltic pump,
- Gear pump,
- Air pressure,
- Dedicated myocardial generator,

- Branching aorta.

Supply

The supply can be realised by a direct connection to the tap (water mains) or via a reservoir which can be filled directly from the tap or manually, e.g. using watering pots. In case of a closed circuit, see section 3.2.5, the reservoir can be filled by the outflow of the phantom.

The D-SPECT room only has a regular indoor faucet making it more difficult to reliably connect to, and the faucet need to be constantly opened and closed between experiments to prevent overflowing. Therefore, although more manual labour, it is safer to manually fill an input reservoir using, for example, watering pots.

Disposal

Proper disposal of contaminated water (nuclear tracer) will be a vital. Any fluid, or materials, that come into contact with the radioactive tracer, needs to be isolated for a period of multiple days. Spills should be avoided at all costs. Waste fluid can be stored in reservoirs, wheeled (closed) containers, or wheelie bins (Dutch "Kliko"). In case of a closed circuit, see section 3.2.5, the outgoing flow can be directed to the input reservoir.

In any case, the contaminated fluid must be collected in a closed container due to safety reasons. The size of the closed container will depend on the configuration and on the physical limitations at the ZGT.

Configuration

The flow set-up can be designed in either a closed circuit or open circuit. The closed circuit characterises itself by having no continuous disposal flow, see figure 3.2. The perfusate from the reservoir must be renewed between experiments. The optional filter can extract, if possible, the tracer from the perfusate such that first pass perfusion is realised. Otherwise, the tracer is recirculated which causes the bolus to disappear. An alternative is an open circuit, see figure 3.3, which guarantees first pass perfusion due to the absence of any form of recirculation.




Figure 3.2: Closed circuit schematic design

Figure 3.3: Open circuit schematic design

Type

As described in section 3.1.4, the flow type is fixed.

Perfusate

Two types of perfusate can be chosen from: water and blood-mimicking fluid (BMF). Water is most practical since it is available in a steady supply. However, BMF is of the same viscosity as human blood allowing for better simulations. However, it needs to be made before any experiments can be performed and in such quantity that it will not run out. In an open circuit

configuration, see section 3.2.5, BMF will rapidly be put into the waste collection making it very costly and would therefore be more suitable in a closed circuit configuration.

3.3 Concepts

Table 3.1: Concept overview

		Concept 1	Concept 2	Concept 3
Myocardium	Shape	VLA: horseshoe, HLA: horseshoe, SA: circle		
	Model	Simplified	Simplified	Simplified
Modality		D-SPECT, dynamic scanning		
Design		LV	LV/RV	LV/RV
Tracer		^{99m}Tc Tetrofosmin via infusion pump		
	Generator	• Gear pumps, • Dmg*.	• Air pressure, • Dmg*.	• Gear pumps, • Branching aorta.
	Supply	Reservoir		
Flow	Disposal	Closed container		
	Configuration	Open	Open	Closed
	Type	Non-pulsatile		
	Perfusate	Water	Water	BMF

* Dedicated myocardial generator

		Concept 4
Myocardium	Shape	VLA: horseshoe, HLA: horseshoe, SA: circle
	Model	17-segments
Modality		D-SPECT, dynamic scanning
Design		4-chamber
Tracer		^{99m}Tc Tetrofosmin via infusion pump
	Generator	• Peristaltic pumps, • Branching aorta.
	Supply	Reservoir
Flow	Disposal	Closed container
	Configuration	Closed
	Type	Pulsatile
	Perfusate	BMF

3.3.1 Concept 1

The characteristic feature of the first concept is its 1-chamber design; only the left ventricle. The ventricle will be surrounded by a fewer-segment myocardium. The simplified nature will make the fabrication more straight-forward. Its ventricle and myocardium are separately supplied using gear pumps for simpler, and more accurate, control of the non-pulsatile flow. The open configuration does not require any filtering but will produce more waste. Gear pumps physically interact with the fluid, allowing for more direct control, which makes water a more suitable choice. Furthermore, it is more practical and less expensive in an open configuration.

216

Pros

- **Simplified myocardium:**
simpler manufacturing.
design+, *fabrication+*
- **1-chamber design:**
simpler manufacturing.
design+, *fabrication+*
- **Straight-forward control:**
more reliable and accurate.
fabrication+
- **Fast response:**
geared pumps respond more direct to control signals making it more flexible if flow resistance is more dynamic.
practicality+
- **Use of water:**
more practical and less expensive.
practicality+
- **Open configuration:**
no filtering.
fabrication+.

Cons

- **Simplified myocardium:**
not completely clinically compatible.
versatility-, *compatibility-*
- **1-chamber design:**
no cardiac noise from surrounding chambers and/or myocardium.
versatility-
- **Ratio between aorta and myocardium flow is manually controlled:**
more difficult to realise physiological effects, i.e. proper contrast ratio.
practicality-, *compatibility-*
- **Open configuration:**
more waste.
practicability-
- **Use of water:**
not a physiological simulation of blood
compatibility-

217 **3.3.2 Concept 2**

218 The characteristic feature of the second concept is its 2-chamber design. Similar to the first
 219 concept, the simplified myocardium allows for simpler manufacturing. The aortic flow flows
 220 through the right ventricle into the left ventricle. The added benefit of the right ventricle, is
 221 that it opens the door to create image noise, i.e. tracer moving through right ventricle or tracer
 222 moving through the right myocardial chamber (if present). In MPI scans, the right myocardium
 223 can clearly be seen, see figure 3.1. As for the generation of the flow, air pressure is used. Flow
 224 remains relatively stable if the pressure remains constant which is both technically and prac-
 225 tically difficult to achieve and less flexible during dynamic flow variations. A compressor could
 226 be a solution but tend to be very bulky and loud, which is not practical in a hospital setting.

227

Pros

- **Simplified myocardium:**
simpler manufacturing.
design+, *fabrication+*
- **2-chamber design:**
ability to create image noise, room for lung-like element.
versatility+, *compatibility+*
- **Air pressure as generator:**
no electromagnetic interference.
other+
- **Use of water:**
more practical and less expensive.
practicality+
- **Open configuration:**
no filtering.
fabrication+

Cons

- **Simplified myocardium:**
not completely clinically compatible.
versatility-, *compatibility-*
- **2-chamber design:**
more difficult to manufacture.
design-, *fabrication-*
- **Air pressure as generator:**
more difficult to control, technically and practically difficult to achieve, not flexible.
fabrication-, *versatility-*
- **Ratio between aorta and myocardium flow is manually controlled:**
more difficult to realise physiological effects, i.e. ratio of tracer that enters the myocardium.
practicality-, *compatibility-*
- **Open configuration:**
more waste.
practicality-
- **Use of water:**
not a physiological simulation of blood
compatibility-

228

229 **3.3.3 Concept 3**

230 The third concept characterises itself by the closed loop configuration. Unlike the first two
 231 concepts, the tracer is removed from the fluid before circulating. It significantly reduces the
 232 amount of waste produced. This enables the use of BMF and the branching aorta. The myocar-
 233 dium is then supplied via a side branch of the aortic flow. The flow to the myocardium can be
 234 regulated by increasing or decreasing the flow resistance in the side branch relative to the aorta.
 235 This technique is theoretically feasible but must be proven. The BMF will more accurately sim-
 236 ulate the behaviour of blood in the phantom but it is costly to make and possibly difficult to
 237 extract the tracer before recirculating. Similar to concept 2, a 2-chamber design is used.

238

Pros

- **Simplified myocardium:**
simpler manufacturing.
design+, *fabrication+*
- **2-chamber design:**
ability to create image noise, room for lung-like element.
versatility+, *compatibility+*
- **Straight-forward control:**
more reliable and accurate.
fabrication+
- **Fast response:**
geared pumps respond more direct to control signals making it more flexible if flow resistance is more dynamic.
practicality+
- **Branched aorta:**
flow ratio is indirectly controlled which makes it easier to realise physiological effects, i.e. ratio of tracer that enters the myocardium.
practicality+, *compatibility+*
- **Closed configuration:**
significant waste reduction.
practicality+
- **Use of BMF:**
more physiological simulation of blood.
compatibility+

239

Cons

- **Simplified myocardium:**
not completely clinically compatible.
versatility-, *compatibility-*
- **2-chamber design:**
more difficult to manufacture.
design-, *fabrication-*
- **Branched aorta:**
unproven.
practicality-
- **Closed configuration:**
requires filtering, can be especially difficult when combined with BMF.
design-, *fabrication-*
- **Use of BMF:**
more costly, difficult to make.
fabrication-, *practicality-*

240 **3.3.4 Concept 4**

241 The fourth concept characterises itself by the 4-chamber design, full myocardium model and,
 242 against system requirements, the pulsatile flow. The 4-chambers simulate the flow through-
 243 out the heart, allowing for mixing the tracer within the chambers. The right ventricle and
 244 left atrium can be directly connected or via an additional element that simulates the lungs.
 245 The full myocardium model, i.e. all 17 segments, ensure that the phantom is compatible
 246 with the clinical practise. The pulsatile flow is a direct consequence of the generator, i.e. the
 247 peristaltic pump. The peristaltic pump has as an added benefit that it is not in direct con-

248 tact with the perfusate making it more suitable to generate flow in combination with BMF.

249 **Pros**

- **17-segment model:**
fully compatible with clinical practice and software.
versatility+, *compatibility+* +
- **4-chamber design:**
full simulation of the heart to simulate the mixing of the tracer in the chambers, room for lung-like element.
versatility+, *compatibility+* +
- **Usage of BMF:**
more physiological simulation of blood.
compatibility+
- **Closed configuration:**
significantly reduced waste.
practicality+
- **Branched aorta:**
flow ratio is indirectly controlled which makes it easier to realise physiological effects, i.e. ratio of tracer that enters the myocardium.
practicality+, *compatibility+*
- **Pulsatile flow:**
more accurate simulation of flow behaviour in human body.
compatibility+

Cons

- **4-chamber design:**
more difficult to manufacture.
design-, *fabrication-*
- **17-segment model:**
more difficult to manufacture.
design-, *fabrication-*
- **Branched aorta:**
unproven.
practicality-
- **Closed configuration:**
requires filtering, can be especially difficult when combined with BMF.
design- *fabrication-*
- **Use of BMF:**
more costly, difficult to make.
fabrication-, *practicality-*
- **Pulsatile flow:**
against system requirements, more difficult to control.
design-, *fabrication-*, *practicality-*, *versatility-*

250 3.4 Comparison

251 3.4.1 Categories

252 The concepts are compared on 5 main categories and a residue category:

- 253 • **Design:** expected design complexity.
254 *A more positive rating indicates a less complicated design of the overall system.*
- 255 • **Fabrication:** expected fabrication complexity.
256 *A more positive rating indicates an overall set-up that is easier to build and program.*
- 257 • **Versatility:** expected flexibility of the system.
258 *A more positive rating indicates that more different scenarios can be simulated without*
259 *major remodelling during experiments; i.e. different kinds of stenosis, generate different*
260 *kinds of noise.*
- 261 • **Practicality:** expected ease of use.
262 *A more positive rating indicates that the system is easier to use; i.e. less continuous waste,*
263 *readily available perfusate, easy to achieve tracer ratio aorta and myocardium.*
- 264 • **Compatibility:** expected equivalence to nature and clinical practise.
265 *A more positive rating indicates a better simulation of the human heart and compatibility*
266 *to clinical protocol and software.*
- 267 • **Other:** other factors that influence the concept's rating.

3.4.2 Weighting factors

Each category has a certain weighting factor by which the pros and cons are multiplied. The combined weights total to 10 points.

- **Design:** receives a weight of 1.
The complexity of the design is not of significant impact assuming it remains feasible.
- **Fabrication:** receives a weight of 3.
Phase 1 is the proof-of-concept phase. The proof-of-concept prototype characterises itself by being relatively simple to fabricate phantom such that more time is available for experiments to gather knowledge and insight for the final phase.
- **Versatility:** receives a weight of 1.
Select experiments are performed which can be prepared beforehand therefore a versatile phantom is not necessary (but can be beneficial).
- **Practicality:** receives a weight of 2.
An easier to use system will ensure that more experiments can be performed in certain amount of time which provides more knowledge and insight.
- **Compatibility:** receives a weight of 2.
The phantom must be compatible with nature and clinical practise but simplifying assumptions are allowed to ease the development of the phase 1 prototype.
- **other:**

3.4.3 Overview

Table 3.2 shows the pros and cons of all the concepts. The blue coloured cells show the most notable characteristics of the concepts; concept 1 excels at the ease of fabrication whilst concept 4 excels at the being true-to-nature and clinically compatible.

Table 3.2: First prototype concept comparison

	Weight	Concept 1	Concept 2	Concept 3	Concept 4
Design	1	++	+-	+-	----
Fabrication	3	++++	+-	+-	-----
Versatility	1	--	+-	+-	++
Practicality	2	+-	+-	+-	+-
Compatibility	2	--	+-	+-	+++++
other	1		+		
Total	10	6+	6-	2+	6-

The total score is calculated using 3.1 where \oplus and \ominus are the + and -, respectively, in table 3.2 for each category n .

$$\sum_{n=design}^{other} (\oplus_n - \ominus_n) * Weight_n \quad (3.1)$$

293 Concept 1 scores highest while concept 2 and 4 have a joint last place. Notably is that concept
294 4 primarily scores poorly on design and fabrication, which are less of an issue in phase 2.

295 **3.5 Definitive choice**

296

4 Detailed design

A Appendix: Mind map




Figure A.1: Mind map for concept designs

Bibliography

Niwaz, T. (2015). Nuclear cardiology: braunwald's chapter 16, disease detection, risk stratification, and clinical decision making.